

Compatibilities



Compatibility for iNedit and double microcatheters

Double microcatheter at iNedit

Microcatheter

Double microcatheter at iNedit			Microcatheter																															
Device	Length (cm)	Maximum OD (in)	Echelon 10	Headway Duo	Phenom 17	Prowler 10	NeuroSlider17	Prowler 14	Echelon 14	Excelsior SL-10	Excelsior HT-17	Headway 17	Headway 21	iNdeep 17	Excelsior 1018	NeuroSlider 21 DLC	Phenom 21	Apollo	Marathon	Trevo trak 21	Via 21	iNdeep 21	Prowler Plus	Excelsior HT-27	Velocity	iNdeep 27	Phenom 27	Headway 27	Via 27	Marksman 027	AHS Offset Delivery	NeuroSlider39		
Echelon 10	150	0.028"	0.002"	0.002"																														
Headway Duo	156/167	0.028"	0.002"	0.002"																														
Phenom 17	150	0.029"																																
Prowler 10	150	0.030"																																
NeuroSlider17	155/160/167	0.030"																																
Prowler 14	150/170	0.030"																																
Echelon 14	150	0.031"																																
Excelsior SL-10	150	0.031"																																
Excelsior HT-17	150	0.031"																																
Headway 17	156	0.032"																																
Headway 21	156	0.033"																																
iNdeep 17	165	0.034"																																
Excelsior 1018	150	0.034"																																
NeuroSlider 21	155/160/167	0.034"																																
Phenom 21	150/160	0.034"																																
Apollo	165	0.035"																																
Marathon	165	0.035"																																
Trevo trak 21	150	0.035"																																
Via 21	154	0.036"																																
iNdeep 21	150/165	0.037"																																
Prowler Plus	150	0.037"																																
Excelsior HT-27	150	0.038"																																
Velocity	160	0.039"																																
iNdeep 27	150/165	0.039"																																
Phenom 27	150/160	0.040"																																
NeuroSlider27	155	0.040"																																
Headway 27	156	0.041"																																
Via 27	135	0.042"																																
Marksman 027	105/135/150/160	0.042"																																
AHS Offset	150	0.050"																																
NeuroSlider39	125/135/145	0.053"																																

✗
Not compatible



≥ 0.003" Compatible



0.002" Limited compatibility



<0.002" Not compatible

Numbers illustrate the remaining ID space within iNedit and two microcatheters



